

## DATASHEET

# DUAL CEM

## Dual Resin Dental Cement

Rebuild lasting

Dual polymerization and dual fixation definitive cement for long-lasting foundations.

**Dual cure, photoactivation, and chemical activation.**

Available in colors A1, A2, A3, and translucent.

### ✓ Indicado para:

- Cementation of crowns, bridges, inlays, onlays, overlays, crowns, and fragments.
- Cementation on zirconia, lithium disilicate, metal porcelain, composite resin, and fiberglass.
- Pre-endodontic restorations.
- Cementation of posts.
- Cementation of implant interfaces.

✓ **Dual-barrel syringe**  
1:1 base and catalyst mixing ratio.

✓ **High radiopacity**  
Enhances material visibility in radiographic images.

✓ **High mechanical properties**

✓ **Easy handling**

✓ **No need for a separate adhesive step**

✓ **Bubble-free**

✓ Available in a dual-barrel syringe of 5 g, with 5 tips for automatic and manual mixing, and instructional manual.

