

## DATASHEET

## K-DAM

**Gingival barrier**

Photopolymerizable composite resin liquid for blocking and protecting tissues during relative or absolute isolation.

It is a direct application product, easy to remove, with excellent thixotropy, and does not drip.

**✓ Indicated for:**

- Complete isolation of soft tissues in dental whitening treatments for clinical use.
- Isolation of soft tissues in dental microabrasion treatments.
- Auxiliary isolation with high sealing power in absolute isolation of isolated parts.

**✓ Excellent sealing power.****✓ High light absorption capacity and dimensional control. Non-irritating to tissues.****✓ Available in a 2g syringe with 5 application tips.**

EXCELLENT  
SEALING POWER.

NON-IRRITATING TO  
GINGIVAL TISSUES.