

K-BLEACH 35%



DATASHEET

K-BLEACH 35%

Drop by drop

Maximum whitening effect without waste

Teeth-whitening kit for clinical use based on 35% hydrogen peroxide with heat blocker for vital and non-vital teeth.

- 18** APPLICATIONS
-  SUITABLE FOR VITAL AND NON-VITAL TEETH
-  HEAT BLOCKER, DOES NOT DEHYDRATE.
-  15-MINUTE SESSIONS WITH COLOR CHANGE



- ✓ **Heat blocker**
The thickener contains a heat blocker to prevent overheating of the tooth in case a lamp is used (light is not necessary).
- ✓ **Drop-by-drop mixing system**
Allows precise dosing for the number of teeth to be whitened, thus avoiding material waste.

PRESENTATION:

- 1 bottle of 8.63ml peroxide.
- 1 bottle of 3.72ml thickener.
- 1 bottle of 1.71ml neutralizer.
- 1 2g K-Dam syringe with application tips.
- 1 spatula.
- 1 mixing plate

K-BLEACH

35%

✓ **Safety mechanism**

Gel color change during the clinical session to maintain a stable pH and prevent enamel demineralization. Initially, the gel is purple and changes to green after a few minutes.

✓ **Excellent viscosity**

Safe application without gel running into unwanted areas.

✓ **Guaranteed safety**

Maximum pH control to prevent enamel demineralization.

✓ **Precise dosage**

Gel savings in each treatment thanks to the drop-by-drop mixing system that allows precise dosing without waste.

✓ **Gel application with the help of a spatula or brush.**



KIYOMI