



KIYOMI

SPECIALISTS IN  
COSMETIC DENTISTRY

# 2024 PRODUCT CATALOGUE



# INDEX

WHITENING

03

ADHESIVES

07

RESINS

09

CEMENTS

14

POLYMERIZATION LAMP

16

REMINERALIZATION, DESSENSITIZING  
& REPAIR

18



WHITER TEETH  
WITH GREAT RESULTS

# K-BLEACH 35%

## CLINIC WHITENING

Tooth whitening based on Hydrogen Peroxide 35% for vital and non-vital teeth in clinic.

- With heat blocker.
- Contains potassium nitrate and sodium fluoride to **reduce hypersensitivity**.
- **Prevents waste:** presentation in two bottles, peroxide and thickener, to mix the necessary amount **in each application**.
- Contains a neutralizing solution bottle in case of contact with soft tissue.





HEAT INHIBITOR,  
DOES NOT DEHYDRATE



15-MINUTE SESSIONS  
WITH COLOR CHANGE



18 APPLICATIONS

### EACH KIT CONTAINS

- 1 bottle of hydrogen peroxide.
- 1 bottle of thickener.
- 1 bottle of neutralizing solution of peroxide (Neutralize).
- 1 spatula.
- 1 plate for gel preparation.
- 1 K-Dam syringe.



# K-BLEACH 16% HOME

## HOME WHITENING

Teeth whitening to be used at home of 16% carbamide peroxide with the supervision of a dentist.

- **Great results** in less time.
- **Reduces dental hypersensitivity** before, during and after treatment.
- Specially indicated for **moderate and severe cases**.
- Its high viscosity and syringe design **facilitate drip-less application**.
- **Prolongs the whitening effects** for longer.

### AVAILABLE IN

- Kit of 4 syringes of 3 ml each.
- Kit of 2 syringes + high viscosity K-Desens 5% dessensitizing gel with potassium nitrate to relieve dental hypersensitivity.





# K-BLEACH RM

## MICROABRASIVE

Microabrasive based on **6% hydrochloric acid** and **silicon carbide** to remove stains due to fluorosis, superficial enamel hypoplasia, structural defects that pigment and white spots of inactive cavities.




REVITALIZES YOUR SMILE



MINIMALLY INVASIVE



GRAY COLOUR TO FACILITATE APPLICATION



ENHANCED MDP



BISPHENOL FREE



5 ML BOTTLE

# K-BOND+ UNIVERSAL

## LIGHT-CURABLE ADHESIVE

Universal **light-curable adhesive** system for enamel and dentin with great adhesion power for longer lasting results.

**Without color alterations**, it does not interfere with the aesthetics of the restoration and final cementation.

- Enhanced MDP: double adhesive power.
- Various engraving protocols: self-etching, engraving selective and total engraving
- Reduces postoperative sensitivity.
- BPA-free.
- It does not interfere with the final result.
- Primer + bond in the same bottle.
- Ethanol base.



#### COLORLESS APPEARANCE

**Does not interfere with the look of restorations and cementations**



#### SMART-LOCK

**Smart double closure for more security and less evaporation**





# ASTER FLOW

## COMPOSITE FLOW

Microhybrid resin with **excellent viscosity**, with a **72% filler charge**, for minimally invasive restorations.

**Reduces microfractures** generated by stress during chewing.

- Excellent fluidity
- High mechanical resistance to chewing.
- Low polymerization shrinkage.



**Aster Flow is indicated for:**

- Sealing pits and fissures.
- Less invasive restorations.
- Base for direct and indirect restorations.
- Prepared in tunnel.
- Repair of enamel defects.
- Restoration in deciduous teeth (baby or temporary teeth).
- Classes III and V.



AVAILABLE COLORS

A1, A2, A3 and B0,5 (2 g syringe)

**72%  
charge**



# ASTER

## COMPOSITE

Microhybrid photocurable resin with nanohybrid particles for **direct and indirect restorations** of anterior and posterior teeth.

Average particle size of 0,5 µm.

- Innovative **silanization** process for greater consistency and better handling.
- Simple and easy color system.
- Natural opalescence and fluorescence.
- Excellent **mechanical properties**.
- Greater polishing and shine.

Aster is specially indicated for **anterior and posterior sectors** (class I, II, III, IV and V).



## FLIP-TOP LID

Easy-to-use lid, it **can be opened and closed with one hand**, and with an airtight seal that guarantees the **conservation of the resin**.



## ERGONOMIC

The syringe body design is **anatomical** so that it fits perfectly in the palm of your hand.

## COLOR FIDELITY

The enamel and dentin shades follow the **VITAPAN Classical scale faithfully**.

## AVAILABLE COLORS

**Dentine DA1, DA2 and DA3 (4 g syringe)**  
**Enamel EA1, EA2 and EA3 (4 g syringe)**  
**B0,5 (2 g syringe)**



# K-DAM

## GINGIVAL BARRIER

Photocurable resin for **total or partial** isolation of the gum.

- Excellent thixotropy, does not drip.
- Direct application.
- Easy to remove.

### **K-DAM is indicated for:**

- Total isolation for in-clinic whitening.
- Partial isolation: dam of isolated pieces.



EXCELLENT  
SEALING POWER

DOES NOT IRRITATE  
GINGIVAL TISSUES



AVAILABLE COLOR  
**BLUE (2 g syringe)**

# DUAL CEM

## DUAL-CURING RESIN CEMENT

**Dual-curing definitive cement** for the cementation of crowns, bridges, fiberglass posts and inlays.

Specially designed for **cementation of lithium, metal and zirconium disilicate.**

- Two types of curing: chemical + photocuring.
- Very resistant adhesion.
- Handleable: contains rheological properties for easy handling.
- No bubbles: the paste mixture does not generate bubbles.
- High mechanical properties.
- Highly radiopaque.



**Dual Cem is indicated for:**

- Cementation of crowns, bridges and inlays.
- Inlays and onlays.
- Laminated crowns and veneers of porcelain or composite resin of hybrid materials.
- Fiberglass, ceramic or metallic posts.
- Maryland Bridges.



WITH EASY  
EXCESS REMOVAL  
SYSTEM



AVAILABLE COLORS

**A1, A2, A3 and TRANS**  
(5 g automix syringe)



**AUTOMIX SYRINGE**

Ensures base mix and catalyst in a 1:1 proportion to increase its adhesive power.



# KIRA

## SMART POLYMERIZATION LAMP

Kira, Kiyomi's polymerization lamp, is the **perfect solution** for your dental practice. It has an **ergonomic and durable design** and a 10 mm diameter metal head, ideal for **light-curing anywhere**, whether posterior or anterior teeth.



### FLAT

Kira is flat, turning her into the perfect tool to work on molars.




### ERGONOMIC DESIGN

Its **compact and ergonomic design** allows it to access any zone of the mouth, in such a way that it ensures that its light beam approaches the photoactivation area as close as possible, achieving a **superior polymerization depth** of the materials.







AMBIDEXTROUS



AUTO POWER OFF

Kira counts with a **smart control chip** to achieve uniform and precise curing. The lamp's 360° head emits light through **4 potent LEDs** (3 blue & 1 ultraviolet) with wave lengths between 385 and 515 nm, with ensures a **wide coverage** for your patients.



ULTRA LIGHTWEIGHT

Weights only 109g.


### Comfort and versatility

The Kira lamp is an **excellent option** for those that seek the **comfort** of a wireless device. The battery **lasts long**, with more than enough energy for your procedures.



Download the KIRA polymerization lamp **catalog** here.





## A HIGHLY EFFECTIVE PREVENTIVE SYMBIOSIS OF HYDROXYAPATITE AND FLUORIDE

### APACARE® REMINERALIZES AND REPAIRS ENAMEL DEFECTS

Liquid enamel, also known as medical hydroxyapatite, contained in ApaCare®, stimulates tooth remineralization, reduces sensitivity and hypersensitivity, and acts as a natural whitening agent.

#### APACARE® REMINERALIZING TOOTHPASTE

SMOOTHENS, SEALS, AND WHITENS TEETH

ApaCare® toothpaste contains sodium fluoride (1450 ppm) and liquid enamel with 3% hydroxyapatite that adheres to the tooth's surface while brushing.



FOR ADULTS AND CHILDREN



PREVENTS CAVITIES

Especially suitable for those prone to cavities, sensitive teeth, or after acidic meals.



### APACARE® LIQUID

SMOOTHER AND BRIGHTER TEETH

**Antibacterial mouthwash with medical hydroxyapatite for caries and periodontitis prevention.**

Alcohol-free, **gentle on the mouth**, and non-irritating.



REDUCES EARLY CAVITIES AND PERIODONTAL BACTERIA



TWICE A DAY

**Contains essential oils that combat biofilms and pathogenic bacteria. Recommended for use twice a day for best results.**



### APACARE® DESSENSITIZING VARNISH HOME

INTENSIVE REPAIR

Dental varnish, ideal for **use with braces**. Prevents cavities and decalcification (white stains).

FULLY BIOCOMPATIBLE PROTECTIVE LAYER



20%

MEDICAL HYDROXYAPATITE


## APACARE® REPAIR

### INTENSIVE REPAIR

ApaCare Repair **desensitizing and repairing gel** is highly enriched with **10% medical hydroxyapatite**.

The liquid dental enamel (medical hydroxyapatite) protects the teeth with a **very thin protective layer**. Teeth are whitened and coated with a shiny layer that **repels bacteria and plaque**. It significantly reduces sensitivity, repairs early caries, pores, and roughness, and **helps slow the progression of white spots**. ApaCare Repair does not contain abrasive substances (RDA <10).


*To be used at home under medical prescription.\**



SIGNIFICANTLY REDUCES SENSITIVITY



HIGHLY ENRICHED WITH HYDROXYAPATITE



DOES NOT CONTAIN ABRASIVE SUBSTANCES (RDA <10)





**APACARE® PROFESSIONAL**  
 INNOVATIVE CLEANING PARTICLES  
 FOR OPTIMAL RESULTS

Polishing paste based on **medical hydroxyapatite** and **sodium fluoride** for professional dental hygiene.



EFFECTIVELY REMOVES  
 PLAQUE AND STAINING



MULTIPLE TREATMENTS

Especially recommended after  
 periodontal therapy and professional  
 dental hygiene.

**APACARE® BARNIZ DESENSIBILIZANTE**  
 INTENSIVE REPAIR

**Natural bonding varnish** for remineralization in cases of sensitive teeth or incipient cavities.

FULLY BIOCOMPATIBLE  
 PROTECTIVE LAYER



**20%** MEDICAL HYDROXYAPATITE





KIYOMI



[www.kiyomidental.com](http://www.kiyomidental.com)



IPG Dental