



KIYOMI

SPECIALISTS IN  
COSMETIC DENTISTRY

# CATALOGUE 2025



# INDEX

WHITENING

03

ADHESIVES

07

COMPOSITE RESINS

09

CEMENTS

14

POLYMERIZATION LAMP

16

REMINERALIZING, DESENSITIZING  
AND REPAIR

18



WHITER TEETH WITH  
GREATER RESULTS

# K-BLEACH 35% CLINIC

## CLINIC TEETH WHITENING

Teeth whitening based on 35% Hydrogen Peroxide for vital and non-vital teeth in clinic.

- With heat blocker
- Contains potassium nitrate and sodium fluoride to **reduce hypersensitivity**.
- **Avoids waste:** presentation in two bottles, peroxide and thickener, to mix the necessary amount in each application.
- It contains a bottle of neutralizer for use in case of contact with soft tissues.





HEAT BLOCKER,  
NO DEHYDRATION



15 MINUTE SESSIONS  
WITH COLOR CHANGE



18 APPLICATIONS

### EACH KIT CONTAINS

- 1 bottle of hydrogen peroxide.
- 1 bottle of thickener.
- 1 bottle of neutralizing solution of peroxide (Neutralize).
- 1 spatula.
- 1 plate for gen preparation.
- 1 K-Dam syringe.



# K-BLEACH 16% HOME

## TEETH WHITENING FOR HOME USE

Teeth whitening to be used at home of 16% carbamide peroxide with the supervision of a dentist.

- **Greater results** in less time.
- **Reduces dental hypersensitivity** before, during and after treatment.
- Specially indicated for **moderate and grave cases**.
- Its high viscosity and syringe design **ease application** without dripping.
- **Prolongs the whitening effects** for much longer.

### AVAILABLE IN

- 4 syringe kit of 3 ml each.
- 2 syringe kit + high viscosity K-Desens 5% desensitizing gel with potassium nitrate to alleviate dental hypersensitivity.





# K-BLEACH RM

## MICROABRASIVE

Microabrasive based on **6% hydrochloric acid and silicon carbide** to remove stains due to fluorosis, superficial enamel hypoplasia, structural defects that pigment and white spots of inactive cavities.




REVITALIZES YOUR SMILE



MINIMALLY INVASIVE



GREY COLOR TO EASE ITS APPLICATION



ENHANCED MDP



FREE OF BISPHENOL



5ML BOTTLE

# K-BOND+ UNIVERSAL

## LIGHT-CURING ADHESIVE

Universal **light-curing adhesive system** for enamel and dentin with great adhesion power for longer lasting results.

**Without color alterations**, it does not interfere with the aesthetics of the restoration and final cementation.

- Enhanced MDP: double adhesive power.
- Various engraving protocols: self-etching, engraving selective and total engraving.
- Reduces postoperative sensitivity.
- BPA-free.
- It does not interfere with the final result.
- Primer + bond in the same bottle.
- Ethanol-based solvent.



### COLORLESS APPEARANCE

**It does not interfere with the aesthetics of restorations and cementations**



### SMART-LOCK

**Smart double closure for more security and less evaporation**





# ASTER FLOW

## FLUID COMPOSITE RESIN

Microhybrid resin with **excellent viscosity**, with a **72% filler load**, for minimally invasive restorations.

**Reduces microfractures** generated by stress during chewing.

- Excellent fluidity.
- High mechanical resistance to chewing.
- Low polymerization shrinkage.



**Aster Flow is indicated for:**

- Sealing pits and fissures.
- Less invasive restorations.
- Base for direct and indirect restorations.
- Prepared in tunnel.
- Repair of enamel defects.
- Restoration in deciduous teeth (baby or temporary teeth).
- Classes III and V.



AVAILABLE COLORS

A1, A2, A3 and B0,5 (2 g syringe)

**72%**  
load



# ASTER

## COMPOSITE RESIN

Microhybrid photocurable resin with nanohybrid particles for **direct and indirect restorations** of anterior and posterior teeth. Average particle size of 0,5µm.

- Innovative **silanization process** for greater consistency and better handling.
- Simple and easy color system
- natural fluorescence and opalescence.
- Excellent **mechanical properties**.
- Even greater polish and shine.

Aster is especially indicated for the **anterior and posterior sector** (class I, II, III, IV and V).



## FLIP-TOP LID

Easy-to-use lid. It **can be opened and closed with one hand**, and with an airtight seal that guarantees the **conservation of the resin**.



## ERGONOMIC

The syringe body design is **anatomical** so that it fits perfectly in the palm of your hand.

## COLOR FIDELITY

The enamel and dentin shades follow the **VITAPAN Classical scale faithfully**.

## AVAILABLE COLORS

**Dentine DA1, DA2 and DA3 (4 g syringe)**  
**Enamel EA1, EA2 and EA3 (4 g syringe)**  
**B0,5 (2 g syringe)**



# K-DAM

## GINGIVAL BARRIER

Photocurable resin for **total or partial isolation** of the gum.

- Excellent thixotropy. Does not drip.
- Direct application.
- Easy to remove.


### K-DAM is indicated for:

- Total isolation for in-clinic teeth whitening
- Partial isolation: dam of isolated pieces.



EXCELLENT  
SEALING POWER

DOES NOT IRRITATE  
GINGIVAL TISSUES.



AVAILABLE COLOR  
**BLUE (2 g syringe)**

# DUAL CEM

## DUAL-CURING CEMENT

Dual-cure permanent cement for the **luting of crowns, bridges, fiberglass posts, and inlays.**

Specifically designed for the cementation of lithium disilicate, metal, zirconia, fiberglass, composite, and hybrid materials.

- Two types of curing: auto & photo.
- High resistance and sealing safety.
- Handleable: contains rheological properties for easy handling.
- High mechanical properties.
- Highly radiopaque.



**Dual Cem is indicated for:**

- Crowns, bridges, posts and inlays.



WITH EASY EXCESS  
REMOVAL SYSTEM



AVAILABLE COLORS

**A1, A2, A3 and TRANS**  
(5 g automix syringe)



**AUTOMIX SYRINGE**

Ensures the base and catalyzer mixture in a 1:1 proportion to increase its adhesive power.



# KIRA

## SMART POLYMERIZATION LAMP

The Kira polymerization lamp by Kiyomi is the **perfect solution** for your dental practice. It features a durable **ergonomic design** and a **10 mm diameter metal head**, making it ideal for light-curing in any area, whether on **anterior teeth or molars**.



### FLAT

Kira is flat, turning it into the perfect tool to work on molars.




### ERGONOMIC DESIGN

Its **ergonomic and compact design** allows access to any area of the mouth, ensuring that the light beam gets as close as possible to the photoactivation area. This results in **greater polymerization depth** of the materials.







AMBIDEXTROUS



AUTO POWER OFF

Kira features an **intelligent control chip** to achieve uniform and precise curing. The lamp's 360° head emits light through **four powerful LEDs** (three blue and one ultraviolet) with wavelengths ranging from 385 to 515 nm, ensuring **broad coverage** for your patients.



ULTRA LIGHTWEIGHT

Weighs just 109g.


#### Comfort and versatility

The Kira lamp is an **excellent choice** for those seeking the **convenience** of a cordless device. Its battery offers **long-lasting performance and durability**, providing sufficient power for your treatments.



Download the **brochure** for the KIRA polymerization lamp here.





## HIGHLY EFFECTIVE PREVENTIVE SYMBIOSIS OF HYDROXYAPATITE AND FLUORIDE

### APACARE® REMINERALIZES AND REPAIRS ENAMEL DEFECT

Liquid enamel, also known as medical hydroxyapatite, which ApaCare® contains, stimulates the remineralization of teeth, reduces sensitivity and hypersensitivity and acts as a natural whitening agent.

#### APACARE® REMINERALIZING PASTE

SMOOTHES, SEALS AND WHITENS TEETH

The ApaCare® toothpaste contains sodium fluoride (1,450 ppm) and liquid enamel with 10% hydroxyapatite, which adheres to the surface of the teeth during brushing.



SUITABLE FOR ADULTS AND CHILDREN



PREVENTS CAVITIES

Especially suitable for individuals with a high tendency to develop cavities, sensitive teeth, or after consuming acidic foods.





## APACARE® LIQUID

SOFTER AND BRIGHTER TEETH

Antibacterial mouthwash with **medical-grade hydroxyapatite** for the **prevention of cavities and periodontitis**.

It is alcohol-free, making it **gentle on the mouth** and preventing irritation.



REDUCES INCIPIENT CAVITIES AND PERIODONTAL BACTERIA



TWICE A DAY

**Contains essential oils that combat biofilms and pathogenic bacteria.**  
It is recommended to use it twice a day for the best results.



## APACARE® VARNISH DESENSITIZING HOME INTENSIVE REPAIR

Dental varnish, ideal **for use with braces**, prevents cavities and demineralization (white spots).

FULLY PROTECTIVE AND BIOCOMPATIBLE LAYER



**20%** MEDICAL HYDROXIAPATITE

## APACARE® REPAIR


### INTENSIVE REPAIR

The ApaCare Repair **restorative and desensitizing gel** is highly enriched with **10% medical hydroxyapatite**.

The liquid enamel (medical hydroxyapatite) protects the teeth with a **very thin protective layer**. It whitens teeth and coats them with a glossy layer that **repels bacteria and plaque**. It significantly reduces sensitivity, repairs early cavities, pores, and rough spots, and **helps slow down the progression of white spots**.

ApaCare Repair does not contain abrasive substances (RDA <10).


*For home use under medical prescription.\**



REDUCES SENSITIVITY  
CONSIDERABLY



HIGHLY ENRICHED WITH  
HIDROXYAPATITE



DOES NOT CONTAIN ABRASIVE  
SUBSTANCES (RDA<10)



### APACARE® PROFESSIONAL

INNOVATIVE CLEANING PARTICLES WITH GREATER RESULTS

Polishing paste based on **medical hydroxyapatite** and **sodium fluoride** for professional dental hygiene.



EFFECTIVELY REMOVES PLAQUE AND STAINING

MULTIPLE TREATMENTS

Especially recommended after periodontal therapy and professional dental hygiene



### APACARE® DESENSITIZING VARNISH

INTENSIVE REPAIR

**Natural adhesive lacquer** for remineralization in case of sensitive teeth or incipient cavities.

FULLY BIOCOMPATIBLE PROTECTIVE LAYER



**20%** MEDICAL HYDROXYAPATITE





KIYOMI



[www.kiyomidental.com](http://www.kiyomidental.com)



IPG Dental