





### **2025 OFFERS**

Valid until December 31, 2025 or end of stock

K-BLEACH 16% HOME

**TEETH WHITENING** FOR HOME USE

K-BLEACH 16% HOME

BRIGHTEN UP YOUR PATIENTS' SMILE



**14**′<sup>53</sup> €

Price of box with 2 g syringe

13<sup>'81</sup> €

Ref. KYMAFLOWBOJ (**B0,5**) Ref. KYMAFLOWAIJ (**A1**) Ref. KYMAFLOWA2J (**A2**) Ref. KYMAFLOWA3J (**A3**)

Unit price with offer

10,36 €

## ASTER FLOW

#### **FLUID COMPOSITE RESIN**

Microhybrid resin with excellent viscosity, with a 72% filler load, for minimally invasive restorations. Reduces microfractures generated by stress during chewing.

- Excellent fluidity.
- High mechanical resistance to chewing..
- Low polymerization shrinkage.

#### **Aster Flow is indicated for:**

- Sealing pits and fissures..
- Less invasive restorations.
- Base for direct restorations.
- Prepared in tunnel.
- Repair of enamel defects.
- Restoration in deciduous tooth (baby or temporary tooth).
- Classes III and V.



Source: Dr. Joaquín Prario



# ASTER

#### NANOHYBRID COMPOSITE RESIN

icrohybrid photocurable resin with nanohybrid particles for direct and indirect restorations of anterior and posterior teeth. Average particle size of 0,5  $\,\mu m$ .

- Innovative silanization process for greater consistency and better handling.
- 4 levels of translucency.
- Opalescence and fluorescence.
- Excellent mechanical properties.
- Polishing and even greater shine.

Aster is especially indicated for the anterior and posterior sector (class I, II, III, IV and V).



**COLOR FIDELITY** 

The enamel and dentin shades follow the VITAPAN Classical scale faithfully.

**19**<sup>'87</sup> €
Price of box
with 4 g syringe

18<sup>'87</sup>

Ref. KYMASDAIJER (DA1) Ref. KYMASDA2JER (DA2) Ref. KYMASDA3JER (DA3) Ref. KYMASEAJJER (EA1) Ref. KYMASEAJJER (EA2) Ref. KYMASEAJJER (EA3)

Unit price with offer

ŧ

**14**'<sup>53</sup> € Price of box with 2 g syringe

13'<sup>81</sup> €

Ref. KYMASB05JER (**B0,5**)


Unit price with offer

10,36 €

**AVAILABLE COLORS:** 

Dentine DA1, DA2 and DA3 (4 g syringe) Enamel EA1, EA2 and EA3 (4 g syringe) B0,5 (2 g syringe)

Source: Dr. Ana Gabriela Alvarado Ramírez



#### FLIP-TOP LID

14,15 €

Easy-to-use lid, it can be opened and closed with one hand, and with an airtight seal that guarantees the conservation of the resin.



#### **ERGONOMIC**

The syringe body design is anatomical so that it fits perfectly in the palm of your hand.



**38**′<sup>37</sup> € Price of box with 5 g automix syringe

**36**′45

Ref. KYMDUALCTRJ (**Trans**) Ref. KYMDUALCAIJ (**A1**) Ref. KYMDUALCA2J (**A2**) Ref. KYMDUALCA3J (**A3**)

Unit price with offer

27,34 €

### WITH EASY EXCESS REMOVAL SYSTEM.

# DUAL

AVAILABLE COLORS:

A1, A2, A3 and TRANS (5 g automix syringe)

#### **DUAL-CURING CEMENT**

**Dual-cure permanent cement** for the luting of crowns, bridges, fiberglass posts, and inlays. Specifically designed for the **cementation of lithium disilicate**, **metal**, **zirconia**, **fiberglass**, **composite**, **and hybrid materials**.

- Two types of curomg: auto & photo.
- High resistence and sealing safety.
- Handleable: contains rheological properties for easy handling.
- High mechanical properties.
- Highly radiopaque.

#### **Dual Cem is indicated for:**

Crowns, bridges, posts and inlays.



AUTOMIX SYRINGE 1:1 proportion

Source: Dr. Leandro Gomila





## K-BOND+ UNIVERSAL

#### **SMART-LOCK**

Smart double closure for more security and less evaporation.

#### LIGHT-CURING ADHESIVE

Universal **light-curing adhesive system** for enamel and dentin with great adhesion power for longer lasting results.

**Without color alterations**, it does not interfere with the aesthetics of the restoration and final cementation.

- Potentialized MDP: double adhesive power.
- Various engraving protocols: self-etching, engraving selective and total engraving.
- Reduces postoperative sensitivity.
- BPA-free.
- It does not interfere with the final result.
- Primer + bond in the same bottle.
- Ethanol-based solvent.



Free of bisphenol



5 ml bottle



Enhanced

MD

COLORLESS APPEARANCE

It does not interfere with the aesthetics of restorations and cementations.



K-BOND+ Universal

Sistema adhesivo autograbante para esmalte y dentina Self-etching adhesive system for enamel and dentin

Adhesivo incoloro | Colorless adhesive

MDP potencializado | Enhanced MPD

Autograbante | Self-etching

**38**′63 €

Price of box with 5 ml bottle

**36**<sup>70</sup>

Ref. KYMKBONDJER

Unit price with offer

27,53 €

Autograbante | Self-d

Source: Dr. Joaquín Prario



#### WHITER TEETH WITH GREATER RESULTS



# K-BLEACH 35% CLINIC



#### **CLINIC TEETH WHITENING**

Tooth whitening based on **35% Hydrogen Peroxide** for vital and non-vital teeth in clinic.

- Comes with a heat blocker.
- Contains potassium nitrate and sodium fluoride to reduce hypersensitivity.
- Avoids waste: presentation in two bottles, peroxide and thickener, to mix the necessary amount in each application.

#### FACH KIT CONTAINS"

- 1 bottle of hydrogen peroxide.
- 1 bottle of thickener.
- 1 bottle of neutralizing solution of peroxide (Neutralize).
- 1 spatula.
- 1 plate for gen preparation.
  - 1 K-Dam syringe.

18 APPLICATIONS

88'<sup>57</sup> €
Kit price

Ref. KYMKBL35KIT

Unit price with offer

63,10 €



Source: Dr. Leandro Gomila







#### AVAILABLE COLOR:

**BLUE (2 g syringe)** 

### K-DAM

#### **GINGIVAL BARRIER**

Photocurable resin for total or partial isolation of the gum.

- Excellent thixotropy. Does not drip.
- Direct application.
- Easy to remove.

#### K-DAM is indicated for:

- Total isolation for in-clinic teeth whitening
- Partial isolation: dam of isolated pieces.

## K-BLEACH RM

STAIN REMOVER BY MICROABRASION

#### REVITALIZE YOUR SMILE

Microabrasive based on 6% hydrochloric acid and silicon carbide to remove stains due to fluorosis, superficial enamel hypoplasia, structural defects that pigment and white spots of inactive cavities.

MINIMALLY INVASIVE

GRAY COLOUR FOR EASY APPLICATION

3+1

24′07 €

Price of box with 2 g syringe

**22**′86

Ref. KYMKBLRMJEF

Unit price with offer

17,15 €









KIT BLANQUEAMIENTO DENTAL

4 jeringas de peróxido de carbamida 16%

# K-BLEACH 16% HOME

3+1

Purchase 3 KITS, receive 1 MORE free of charge (the one with the lowest price)

#### 36'33€

Kit price: 4 syringes of 3 ml

**34**′<sup>52</sup>

Ref. KYMKBL16KT4

Unit price with offer

25,90 €

### TEETH WHITENING FOR HOME USE

Teeth whitening to be used at home of 16% carbamide peroxide with the supervision of a dentist.

- Great results in less time.
- Available in a kit of 4 syringes of 3 ml each and in a kit of 2 syringes of 3 ml + 1 syringe of K-Desens 5% of 3 ml.





Kit price: 2 syringes of 3 ml + 1 syringe of K-Desens 5% of 3 m

**31**′03

€

Ref. KYMKBL16KTD

Unit price with offer

23,27 €

# 3+1

#### 16'33 € Price of 3ml syringe

**15**′<sup>52</sup> €

Ref. KYMKFDESEN5

Unite price with offer

11,64 €

### K-DESENS 5%

#### **DESENSITIZING GEL**

Low viscosity gel with potassium nitrate, designed to relieve dental hypersensitivity.

 To treat tooth sensitivity: before, during or after whitening, and in cases of gingival recession or periodontal treatment.







Renowned for its **advanced ergonomics** and notable **durability**, Kiyomi's Kira curing light provides comfortable and highly efficient handling. Its metal head, with an **ideal diameter of 10 mm** is excellently designed to perform photopolymerization in any area of the mouth, whether on front teeth or molars. Kira's **precision and effectiveness** are essential to ensure the highest quality of dental treatments.

# OFFER! Purchase €1,500 worth of KIYOMI products and receive a free Kira lamp

400<sup>20</sup>€
Unit price

380<sup>20</sup>€











ERGONOMIC DESIGN

ULTRA LIGHTWFIGHT

AUTO POWER OFF

AMBIDEXTROUS

LONG-LASTING





KIRA IS FLAT, TURNING HER INTO THE **PERFECT TOOL**TO WORK ON MOLARS.

#### MULTIPLE CURING OPTIONS

Kira offers six modes of operation to adapt to various clinical situations:





#### MODES MO, M1 & M2

Light intensity adjustment between 1,000 mW/cm $^2$ , 1,500 mW/cm $^2$  or 2,300 mW/cm $^2$ , with a wavelength of 385 Nm to 515 Nm, covering a wide spectrum.

#### **RAMP MODE**

It gradually increases the light intensity from 0 to 1,000mW/cm2 during the first 5 seconds and then remains constant by emitting a sound every 5 seconds.





#### **PULSE MODE**

It emits a **light intensity** of 1,000mW/cm2 that flashes **once per second and emits a sound every 5 seconds**.

#### **DETECTION MODE**

Its light intensity is 600mW/cm2 with a wavelength between 385nm and 410nm. When you press the activation button, the light activates immediately and emits a sound every 5 seconds.



#### **TECHNICAL SPECS**

**Size:** 24x21x230 mm

Net weight: 109 g

Lens size: 11 mm

Power supply:

1200 times, 3 sec maximum light





### HIGHLY EFFECTIVE PREVENTIVE SYMBIOSIS OF HYDROXYAPATITE AND FLUORIDE

#### APACARE® REMINERALIZES AND REPAIRS ENAMEL DEFECTS

Liquid enamel, also known as medical hydroxyapatite, which ApaCare® containins, stimulates the remineralization of teeth, reduces sensitivity and hypersensitivity and acts as a natural whitening agent.

#### **APACARE® REMINERALIZING PASTE**

### SMOOTHES, SEALS AND WHITENS TEETH

ApaCare® toothpaste contains sodium fluoride (1,450 ppm fluorine ion) and liquid enamel that adhere to the surface of the teeth during brushing.



SUITABLE FOR ADULTS AND CHILDREN.



CONTAINS ESSENTIAL OILS THAT FIGHT BIOFILMS AND PATHOGENIC BACTERIA..

### **APACARE® VARNISH DESENSITIZING HOME** INTENSIVE



RFPAIR

**FULLY PROTECTIVE** AND BIOCOMPATIBLE LAYER.



- Application brush integrated into the cap
- for easy use.

7'85**€**. 75 ml tube.

Ref. 1001619



Ref. 1040200







CONTAINS ESSENTIAL OILS THAT FIGHT BIOFILMS AND PATHOGENIC

BACTERIA.

12′67€ 200 ml bottle.



Ref. 10011402



CONSIDERABLY REDUCES SENSITIVITY.



HIGHLY ENRICHED WITH 10% HYDROXYAPATITE.



DOES NOT CONTAIN ABRASIVE SUBSTANCES (RDA<10).



# APACARE® REPAIR INTENSIVE REPAIR

#### Repairing, desensitizing and whitening gel.

ApaCare Repair is highly enriched with medical hydroxyapatite (10%).

The liquid dental enamel (medical hydroxyapatite) protects teeth with an ultra-thin protective layer. It whitens teeth and coats them with a shiny layer that repels bacteria and plaque. It significantly reduces sensitivity, repairs early-stage cavities, pores, and roughness, and helps slow the progression of white spots.

ApaCare Repair contains no abrasive substances (RDA <10).

Home use under medical prescription. \*



Ref. 1001620

15<sup>73</sup>



#### **APACARE® PROFESSIONAL**

### INNOVATIVE CLEANING PARTICLES FOR BIG RESULTS

Polishing paste based on medical hydroxyapatite and sodium fluoride for professional dental hygiene.



EFFECTIVELY REMOVES PLAQUE AND STAINING





ESPECIALLY
RECOMMENDED
AFTER PERIODONTAL
THERAPY AND
PROFESSIONAL
DENTAL HYGIENE.

**36**′85 € 75 ml tube.

**33**′¹7





### INTENSIVE REPAIR

Natural adhesive lacquer for remineralization in case of sensitive teeth or incipient cavities..



**23**′89 € 5 ml bottle.

**21**′50





FULLY BIOCOMPATIBLE PROTECTIVE LAYER.



MEDICAL HYDROXYAPATITE.

